

MOJOFEET

made to move

HOW TO FIT YOUR PATIENT'S MOJOFEET PREFABS.

1. Remove the insole from their current shoes and place it on the bottom of the MojoFeet.
2. Trace the insole onto the MojoFeet.
3. Cut the MojoFeet to size.
4. You may have to trim and shorten the MojoFeet to get a perfect fit.



DON'T BE FOOLED BY SIZES.



The sizing of the MojoFeet is actually best figured out by matching the prefab up to the patient's foot. The big toe should sit where the arch contour ends.

We call this the 'heel to ball of foot size' and it can vary greatly from what the patient actually tells you is their shoe size. Patients should be wearing shoes that flex at this same point (where the big toe flexes). Many times patients are in the wrong shoe size!

HOW TO MOLD YOUR PATIENT'S MOJOFEET PREFABS.

To get a 'more comfortable fit' for patients:

1. WEAR THEM – after a few days they will mold to their feet.
2. HEAT THEM – for a higher level of comfort quicker. Preheat the oven to 200F/90C. (yes, it is time for a toaster oven in your practice!) DO NOT MICROWAVE. Put your MojoFeet on a baking sheet and bake for 2 minutes. After 2 minutes, put them back in your footwear and stand on them right away. Make sure your heels are snug against the back of your shoes and lace them up; standing on your warm and toasty MojoFeet Prefabs for 2 minutes. MojoFeet Prefabs are dispensed in a medical setting, please heat your patient's devices and do not ask them to do it at home.

